



Groynne Tech
be in vision

Specialized in Design of Containment System



GT-RTPLT - Rapid Transfer Leak Test

The GT-RTPLT allows the leak test in less than 5 minutes of a RTP container. The test can be performed prior the container connection to the isolator without breaking the containers integrity.

Technical Details:

Allows test of RTP containers of every diameter from 105 to 350mm

Leak test performed with high pressure decay method
Test records with operator and container identification
test data and result.

Optional Wireless Technology



“Innovation in Technology” is the Motto of the Organisation